




Smarter Thinking. Simpler Design.



Genesis

Product Catalog

Contact Us:

RUNDAS GmbH
Amalienstr. 62
46537 Dinslaken
Tel.: 02064 625 95 50
Fax: 02064 625 95 80
Email: info@rundas.de
Internet: www.rundas.de



TABLE OF CONTENTS

Genesis

Ø3.5/3.8 mm Genesis Implants	4
Ø3.5/3.8 mm Prosthetic Options	5
Ø3.5/3.8 mm Prosthetic Components	6
Ø4.5 mm Genesis Implants	10
Ø4.5 mm Prosthetic Options	11
Ø4.5 mm Prosthetic Components	12
Ø5.5/6.5 mm Genesis Implants	16
Ø5.5/6.5 mm Prosthetic Options	17
Ø5.5/6.5 mm Prosthetic Components	18

Digital

Imagine™ Milling Services	23
Core3dcentres®	26
In-Office Milling	28

LOCATOR® Abutments

LOCATOR® Prosthetic Options	30
LOCATOR® Abutment	31
LOCATOR® R-Tx™ Abutment	32
LOCATOR® F-Tx™ Abutment	33



Surgical Kit & Instruments

Genesis Surgical Kit	35
Genesis Instruments	36


Product specifications are subject to change without notice.
Items illustrated are not to scale.
Certain products may not be available in all countries.



Genesis Implants

Internal TiLobe®

Genesis Implants



Catalog

G21050
G21051
G21052
G21053

G21130
G21131
G21132
G21133
G21135

Implant Size

Ø3.5 x 10 mm Length
Ø3.5 x 11.5 mm Length
Ø3.5 x 13 mm Length
Ø3.5 x 16 mm Length

Ø3.8 x 8.5 mm Length
Ø3.8 x 10 mm Length
Ø3.8 x 11.5 mm Length
Ø3.8 x 13 mm Length
Ø3.8 x 16 mm Length

Cover Screw



Catalog

G21100

Description

Ø3.5/3.8 mm

Healing Abutments



Catalog

G21103
G21104
G21105
G21106

Description

Ø3.5/3.8 mm
Ø3.5/3.8 mm
Ø3.5/3.8 mm
Ø3.5/3.8 mm


Flare

4.0 mm
4.0 mm
5.0 mm
5.0 mm

Cuff Height

3.0 mm
5.0 mm
3.0 mm
5.0 mm

The Ø3.5 & Ø3.8 mm Genesis can be restored with Ø4.0 and Ø5.0 mm flared prosthetic components. An Abutment Screw is included with all TiLobe® abutments.





Prosthetic Components

Internal TiLobe®

Open & Closed Tray Impression Post



Catalog #

Description

Flare

G71100

Ø3.5/3.8 mm

4.0 mm

G71101

Ø3.5/3.8 mm

5.0 mm

Open Tray Impression Post



Catalog #

Description

Flare

45904K

Ø3.5/3.8 mm

4.0 mm

45906K

Ø3.5/3.8 mm

5.0 mm

Closed Tray Impression Post



Catalog #

Description

Flare

45905K

Ø3.5/3.8 mm

4.0 mm

45907K

Ø3.5/3.8 mm

5.0 mm

Implant Analog



Catalog #

Description

G71107

Ø3.5/3.8 mm, 1-Pack

G71108

Ø3.5/3.8 mm, 10-Pack

Temporary Abutment



Catalog #

Description

Flare

G71130

Ø3.5/3.8 mm, Locking

4.0 mm

G71131

Ø3.5/3.8 mm, Non-Locking

4.0 mm

G71132

Ø3.5/3.8 mm, Locking

5.0 mm

G71133

Ø3.5/3.8 mm, Non-Locking

5.0 mm

Ti Temporary Abutment



Catalog #

Description

G71144

Ø3.5/3.8 mm, Locking

G71145

Ø3.5/3.8 mm, Non-Locking

Esthetic Contour Ti Abutment



Catalog #

Description

Flare

Cuff Height

G71200

Ø3.5/3.8 mm, Straight

4.0 mm

1.0 mm

G71201

Ø3.5/3.8 mm, Straight

4.0 mm

2.0 mm

G71202

Ø3.5/3.8 mm, 15° Angled

4.0 mm

1.0 mm

G71203

Ø3.5/3.8 mm, 15° Angled

4.0 mm

2.0 mm

G71204

Ø3.5/3.8 mm, Straight

5.0 mm

1.0 mm

G71205

Ø3.5/3.8 mm, Straight

5.0 mm

3.0 mm

G71206

Ø3.8 mm, 15° Angled

5.0 mm

1.0 mm

G71207

Ø3.8 mm, 15° Angled

5.0 mm

3.0 mm

UCLA Abutment



Catalog #

Description

Cuff Height

G71231

Gold/Plastic Sleeve, Locking

1.0 mm

G71232

Gold/Plastic Sleeve, Non-Locking

1.0 mm

Replacement
Final Abutment Screw



Catalog #

Description

G71125

1-Pack

Laboratory Screw



Catalog #

Description

G71127

10-Pack

Quick-Packs



Catalog #

Description

Flare

Cuff Height

G71400

Ø3.5/3.8 mm

5.0 mm

1.0 mm

G71401

Ø3.5/3.8 mm

5.0 mm

2.0 mm

G71402

Ø3.5/3.8 mm

5.0 mm

3.0 mm


G71403

Ø3.5/3.8 mm

5.0 mm

4.0 mm

Quick-Abutment



Catalog #

Description

Flare

Cuff Height

G71434

Ø3.5/3.8 mm

5.0 mm

1.0 mm

G71435

Ø3.5/3.8 mm

5.0 mm

2.0 mm

G71436

Ø3.5/3.8 mm

5.0 mm

3.0 mm

G71437

Ø3.5/3.8 mm

5.0 mm

4.0 mm

Quick-Abutment
Quick-Cap®



Catalog #

Description

Flare

G71428

Ø3.5/3.8 mm, 1-Pack

5.0 mm

Quick-Abutment
Temporary/Healing Cap



Catalog #

Description

Flare

G71462

Ø3.5/3.8 mm

5.0 mm

Quick-Abutment Analog



Catalog #

Description

Flare

G71468

Ø3.5/3.8 mm

5.0 mm

Quick-Abutment
Plastic Sleeve



Catalog #

Description

Flare

G71478

Ø3.5/3.8 mm, Locking

5.0 mm

G71479

Ø3.5/3.8 mm, Non-Locking

5.0 mm

Quick-Abutment
Ceramic Coping



Catalog #

Description

Flare

G71475

Ø3.5/3.8 mm

5.0 mm



Prosthetic Components

Internal TiLobe®

Multi Unit Abutment



Catalog #

45312K
45313K
45314K

Description

Ø3.5/3.8 mm, Straight
Ø3.5/3.8 mm, Straight
Ø3.5/3.8 mm, Straight

Cuff Height

1.0 mm
2.0 mm
3.0 mm

45321K
45322K

Ø3.5/3.8 mm, Angled, 17°
Ø3.5/3.8 mm, Angled, 17°

3.0 mm
4.0 mm

45936K
45937K
45938K

Ø3.5/3.8 mm, Angled, 30°
Ø3.5/3.8 mm, Angled, 30°
Ø3.5/3.8 mm, Angled, 30°

3.0 mm
4.0 mm
5.0 mm

Multi Unit Abutment
Healing Cap



Catalog #

45330K

Description

Ø3.5/3.8 mm

Multi Unit Abutment
Impression Post



Catalog #

45327K
45328K

Description

Open Tray, Ø3.5/3.8 mm with screw
Closed Tray, Ø3.5/3.8 mm

Multi Unit Abutment
Analog



Catalog #

45329K

Description

Ø3.5/3.8 mm



Multi Unit Abutment
Ti Temporary Sleeve



Catalog #

45331K

Description

Ø3.5/3.8 mm

Multi Unit Abutment
Plastic Sleeve



Catalog #

45332K

Description

Ø3.5/3.8 mm

Multi Unit Abutment
Prosthetic Screw



Catalog #

45333K

Description

1-Pack

Angled Multi-Unit Abutment
Replacement Screw



Catalog #

45334K

Description

1-Pack



Genesis Implants
Internal TiLobe®

Genesis Implants



Catalog #

G21137
G21138
G21139
G21140
G21142

Implant Size

Ø4.5 x 8.5 mm Length
Ø4.5 x 10 mm Length
Ø4.5 x 11.5 mm Length
Ø4.5 x 13 mm Length
Ø4.5 x 16 mm Length

Cover Screw



Catalog #

G21101

Description

Ø4.5 mm

Healing Abutments



Catalog #

G21107
G21108

G21109
G21110

Description

Ø4.5 mm
Ø4.5 mm

Ø4.5 mm
Ø4.5 mm

Flare

5.0 mm
5.0 mm


6.0 mm
6.0 mm

Cuff Height

3.0 mm
5.0 mm

3.0 mm
5.0 mm

The Ø4.5 mm Genesis can be restored with Ø5.0 and Ø6.0 mm flared prosthetic components.
An Abutment Screw is included with all TiLobe® abutments.



Open & Closed Tray
Impression Post

Catalog #

Description

Flare

G71102
G71103

Ø4.5 mm
Ø4.5 mm

5.0 mm
6.0 mm

Open Tray
Impression Post

Catalog #

Description

Flare

45908K
45910K

Ø4.5 mm
Ø4.5 mm

5.0 mm
6.0 mm

Closed Tray
Impression Post

Catalog #

Description

Flare

45909K
45911K

Ø4.5 mm
Ø4.5 mm

5.0 mm
6.0 mm

Implant Analog



Catalog #

Description

G71109
G71110

Ø4.5 mm, 1-Pack
Ø4.5 mm, 10-Pack

Temporary Abutment



Catalog #

Description

Flare

G71134
G71135
G71136
G71137

Ø4.5 mm, Locking
Ø4.5 mm, Non-Locking
Ø4.5 mm, Locking
Ø4.5 mm, Non-Locking

5.0 mm
5.0 mm
6.0 mm
6.0 mm

Ti Temporary Abutment



Catalog #

Description

G71146
G71147

Ø4.5 mm, Locking
Ø4.5 mm, Non-Locking

Esthetic Contour
Ti Abutment

Cuff Height —
(Buccal)

Catalog #

Description

Flare

Cuff Height

G71208
G71209
G71210
G71211

Ø4.5 mm, Straight
Ø4.5 mm, Straight
Ø4.5 mm, 15° Angled
Ø4.5 mm, 15° Angled

5.0 mm
5.0 mm
5.0 mm
5.0 mm

0.5 mm
1.5 mm
0.5 mm
1.5 mm

G71212
G71213
G71214
G71215

Ø4.5 mm, Straight
Ø4.5 mm, Straight
Ø4.5 mm, 15° Angled
Ø4.5 mm, 15° Angled

6.0 mm
6.0 mm
6.0 mm
6.0 mm

1.0 mm
3.0 mm
1.0 mm
3.0 mm

UCLA Abutment



Catalog #

Description

Cuff Height

G71233
G71234

Gold/Plastic Sleeve, Locking
Gold/Plastic Sleeve, Non-Locking

1.0 mm
1.0 mm

Replacement
Final Abutment Screw



Catalog #

Description

G71125

1-Pack

Laboratory Screw



Catalog #

Description

G71127

10-Pack

Quick-Packs



Catalog #

Description

Flare

Cuff Height

G71404

Ø4.5 mm

5.0 mm

0.5 mm

G71405

Ø4.5 mm

5.0 mm

1.5 mm

G71406

Ø4.5 mm

5.0 mm

2.5 mm

G71407

Ø4.5 mm

5.0 mm

4.0 mm

G71408

Ø4.5 mm

6.0 mm

0.5 mm

G71409

Ø4.5 mm

6.0 mm

1.5 mm

G71410

Ø4.5 mm

6.0 mm

2.5 mm

G71411

Ø4.5 mm

6.0 mm

4.0 mm

Quick-Abutment



Catalog #

Description

Flare

Cuff Height

G71438

Ø4.5 mm

5.0 mm

0.5 mm

G71439

Ø4.5 mm

5.0 mm

1.5 mm

G71440

Ø4.5 mm

5.0 mm

2.5 mm

G71441

Ø4.5 mm

5.0 mm

4.0 mm

G71442

Ø4.5 mm

6.0 mm

0.5 mm

G71443

Ø4.5 mm

6.0 mm

1.5 mm

G71444

Ø4.5 mm

6.0 mm

2.5 mm

G71445

Ø4.5 mm

6.0 mm

4.0 mm

Quick-Abutment
Quick-Cap®



Catalog #

Description

Flare

G71429

Ø4.5 mm, 1-Pack

5.0 mm

G71430

Ø4.5 mm, 1-Pack

6.0 mm

Quick-Abutment
Temporary/Healing Cap



Catalog #

Description

Flare

G71463

Ø4.5 mm

5.0 mm

G71464

Ø4.5 mm

6.0 mm

Quick-Abutment
Analog



Catalog #

Description

Flare

G71469

Ø4.5 mm

5.0 mm

G71470

Ø4.5 mm

6.0 mm

Quick-Abutment
Plastic Sleeve



Catalog #

Description

Flare

G71480

Ø4.5 mm, Locking

5.0 mm

G71481

Ø4.5 mm, Non-Locking

5.0 mm

G71482

Ø4.5 mm, Locking

6.0 mm

G71483

Ø4.5 mm, Non-Locking

6.0 mm

Quick-Abutment
Ceramic Coping



Catalog #

Description


Flare

G71476

Ø4.5 mm

5.0 mm

Multi Unit Abutment



Catalog #

Description

Cuff Height

45315K

Ø4.5 mm, Straight

1.0 mm

45316K

Ø4.5 mm, Straight

2.0 mm

45317K

Ø4.5 mm, Straight

3.0 mm

45318K

Ø4.5 mm, Straight

4.0 mm

45323K

Ø4.5 mm, Angled, 17°

3.0 mm

45324K

Ø4.5 mm, Angled, 17°

4.0 mm

45965K

Ø4.5 mm, Angled, 30°

3.0 mm

45325K

Ø4.5 mm, Angled, 30°

4.0 mm

45326K

Ø4.5 mm, Angled, 30°

5.0 mm

Multi Unit Abutment
Healing Cap



Catalog #

Description

45330K

Ø4.5 mm

Multi Unit Abutment
Impression Post



Catalog #

Description

45327K

Open Tray, Ø4.5 mm with screw

45328K

Closed Tray, Ø4.5 mm

Multi Unit Abutment
Analog



Catalog #

Description

45329K

Ø4.5 mm

Multi Unit Abutment
Ti Temporary Sleeve



Catalog #

45331K

Description

Ø4.5 mm

Multi Unit Abutment
Plastic Sleeve



Catalog #

45332K

Description

Ø4.5 mm

Multi Unit Abutment
Prosthetic Screw



Catalog #

45333K

Description

1-Pack

Angled Multi-Unit Abutment
Replacement Screw



Catalog #


45334K

Description

1-Pack




Genesis Implants Internal TiLobe®



Catalog #	Implant Size
G21144	Ø5.5 x 8.5 mm Length
G21145	Ø5.5 x 10 mm Length
G21146	Ø5.5 x 11.5 mm Length
G21147	Ø5.5 x 13 mm Length

Catalog #	Implant Size
G21150	Ø6.5 x 8.5 mm Length
G21151	Ø6.5 x 10 mm Length
G21152	Ø6.5 x 11.5 mm Length
G21153	Ø6.5 x 13 mm Length

Cover Screw




Catalog #	Description
G21102	Ø5.5/6.5 mm

Healing Abutments



Catalog #	Description	Flare	Cuff Height
G21111	Ø5.5 mm	5.0 mm	3.0 mm
G21112	Ø5.5 mm	5.0 mm	5.0 mm
G21113	Ø5.5/6.5 mm	6.0 mm	3.0 mm
G21114	Ø5.5/6.5 mm	6.0 mm	5.0 mm
G21296	Ø5.5/6.5 mm	7.0 mm	3.0 mm
G21297	Ø5.5/6.5 mm	7.0 mm	5.0 mm
G21298	Ø5.5/6.5 mm	8.0 mm	3.0 mm
G21299	Ø5.5/6.5 mm	8.0 mm	5.0 mm

The Ø5.5 and Ø6.5 mm Genesis can be restored with Ø5.0 and Ø6.0 mm flared prosthetic components. An Abutment Screw is included with all TiLobe® abutments.





Prosthetic Components

Internal TiLobe®

Open & Closed Impression Post



Catalog #

Description

Flare

G71104

Ø5.5 mm

5.0 mm

G71105

Ø5.5/6.5 mm

6.0 mm

Impression Post



Catalog #

Description

Flare

45916K

Ø5.5/6.5 mm

7.0 mm

45918K

Ø5.5/6.5 mm

8.0 mm

Impression Post



Catalog #

Description

Flare

45917K

Ø5.5/6.5 mm

7.0 mm

45919K

Ø5.5/6.5 mm

8.0 mm

Implant Analog



Catalog #

Description

G71111

Ø5.5 mm, 1-Pack

G71112

Ø5.5 mm, 10-Pack

G71113

Ø6.5 mm 1-Pack

G71114

Ø6.5 mm 10-Pack

Temporary Abutment



Catalog #

Description

Flare

G71138

Ø5.5 mm, Locking

5.0 mm

G71139

Ø5.5 mm, Non-Locking

5.0 mm

G71140

Ø5.5/6.5 mm, Locking

6.0 mm

G71141

Ø5.5/6.5 mm, Non-Locking

6.0 mm

Ti Temporary Abutment



Catalog #

Description

G71148

Ø5.5/6.5 mm, Locking

G71149

Ø5.5/6.5 mm, Non-Locking

Esthetic Contour Ti Abutment



Catalog #

Description

Flare

Cuff Height

G71216

Ø5.5 mm, Straight

5.0 mm

0.5 mm

G71217

Ø5.5 mm, Straight

5.0 mm

1.5 mm

G71218

Ø5.5 mm, 15° Angled

5.0 mm

0.5 mm

G71219

Ø5.5 mm, 15° Angled

5.0 mm

1.5 mm

G71220

Ø5.5/6.5 mm, Straight

6.0 mm

1.0 mm

G71221

Ø5.5/6.5 mm, Straight

6.0 mm

3.0 mm

G71222

Ø5.5/6.5 mm, 15° Angled

6.0 mm

1.0 mm

UCLA Abutment



Catalog #

Description

Cuff Height

G71235

Gold/Plastic Sleeve, Locking

1.0 mm

G71236

Gold/Plastic Sleeve, Non-Locking

1.0 mm

Replacement
Final Abutment Screw



Catalog #	Description
G71125	1-Pack

Laboratory Screw




Catalog #	Description
G71127	10-Pack

Quick-Packs



Catalog #	Description	Flare	Cuff Height
G71412	Ø5.5 mm	5.0 mm	0.5 mm
G71413	Ø5.5 mm	5.0 mm	1.5 mm
G71414	Ø5.5 mm	5.0 mm	2.5 mm
G71415	Ø5.5 mm	5.0 mm	4.0 mm
G71416	Ø5.5/6.5 mm	6.0 mm	0.5 mm
G71417	Ø5.5/6.5 mm	6.0 mm	1.5 mm
G71418	Ø5.5/6.5 mm	6.0 mm	2.5 mm
G71419	Ø5.5/6.5 mm	6.0 mm	4.0 mm

Quick-Abutment



Catalog #	Description	Flare	Cuff Height
G71446	Ø5.5 mm	5.0 mm	0.5 mm
G71447	Ø5.5 mm	5.0 mm	1.5 mm
G71448	Ø5.5 mm	5.0 mm	2.5 mm
G71449	Ø5.5 mm	5.0 mm	4.0 mm
G71450	Ø5.5/6.5 mm	6.0 mm	0.5 mm
G71451	Ø5.5/6.5 mm	6.0 mm	1.5 mm
G71452	Ø5.5/6.5 mm	6.0 mm	2.5 mm
G71453	Ø5.5/6.5 mm	6.0 mm	4.0 mm

Quick-Abutment
Quick-Cap®



Catalog #	Description	Flare
G71431	Ø5.5 mm, 1-Pack	5.0 mm
G71432	Ø5.5/6.5 mm, 1-Pack	6.0 mm

Quick-Abutment
Temporary/Healing Cap



Catalog #	Description	Flare
G71465	Ø5.5 mm	5.0 mm
G71466	Ø5.5/6.5 mm	6.0 mm

Quick-Abutment Analog



Catalog #	Description	Flare
G71471	Ø5.5 mm	5.0 mm
G71472	Ø5.5/6.5 mm	6.0 mm



Prosthetic Components

Internal TiLobe®

Quick-Abutment Plastic Sleeve




<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
G71484	Ø5.5 mm, Locking	5.0 mm
G71485	Ø5.5 mm, Non-Locking	5.0 mm
G71486	Ø5.5/6.5 mm, Locking	6.0 mm
G71488	Ø5.5/6.5 mm, Non-Locking	6.0 mm
45939K	Ø5.5/6.5 mm, Straight	3.0 mm
45940K	Ø5.5/6.5 mm, Straight	4.0 mm
45941K	Ø5.5/6.5 mm, Angled, 17°	3.0 mm
45942K	Ø5.5/6.5 mm, Angled, 17°	4.0 mm
45966K	Ø5.5/6.5 mm, Angled, 30°	3.0 mm
45967K	Ø5.5/6.5 mm, Angled, 30°	4.0 mm

Quick-Abutment Ceramic Coping



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
G71477	Ø5.5 mm	6.0 mm

Multi Unit Abutment



<i>Catalog #</i>	<i>Description</i>	<i>Cuff Height</i>
45319K	Ø5.5/6.5 mm, Straight	1.0 mm
45320K	Ø5.5/6.5 mm, Straight	2.0 mm
45939K	Ø5.5/6.5 mm, Straight	3.0 mm
45940K	Ø5.5/6.5 mm, Straight	4.0 mm
45941K	Ø5.5/6.5 mm, Angled, 17°	3.0 mm
45942K	Ø5.5/6.5 mm, Angled, 17°	4.0 mm
45966K	Ø5.5/6.5 mm, Angled, 30°	3.0 mm
45967K	Ø5.5/6.5 mm, Angled, 30°	4.0 mm

Multi Unit Abutment Healing Cap



<i>Catalog #</i>	<i>Description</i>
45330K	Ø5.5/6.5 mm

**Multi Unit Abutment
Impression Post**



Catalog #

Description

45327K
45328K

Open Tray, Ø5.5/6.5 mm with screw
Closed Tray, Ø5.5/6.5 mm

**Multi Unit Abutment
Analog**



Catalog #

Description

45329K

Ø5.5/6.5 mm

**Multi Unit Abutment
Ti Temporary Sleeve**



Catalog #

Description

45331K

Ø5.5/6.5 mm

**Multi Unit Abutment
Plastic Sleeve**



Catalog #

Description

45332K

Ø5.0/5.5/6.5 mm

**Multi Unit Abutment
Prosthetic Screw**



Catalog #

Description

45333K

1-Pack

**Angled Multi-Unit Abutment
Replacement Screw**



Catalog #

Description

45334K

1-Pack

Digital

Angled Abutment Drivers
(Latch Type) (Torque Wrench)



Angulated Screw




Ti-Base



Angulated Screw
Channel

Quad Driver
(Torque Wrench) Quad Driver
(Swivel Head)



Final Abutment Screw



Ti-Base



Locking



Non-Locking

Ti-Blanks



Digital Scanbodies



Digital Analog



Download Digital Libraries from Imagine

3shape

exocad

Temporary Abutments



Acrylic



Titanium



Milling Services

Healing Abutment



Ø4.1 Shown



Digital Components

Internal TiLobe®

Digital Scanbodies



Catalog

Description

IM-SBKEP35

Mist IC Scanbody 3.5

IM-SBKEP41

Mist IC Scanbody 4.1

IM-SBKEP50

Mist IC Scanbody 5.0

Digital Analog



Catalog

Description

IM-3DKEP35

Mist IC 3D Analog 3.5

IM-3DKEP41

Mist IC 3D Analog 4.1

IM-3DKEP50

Mist IC 3D Analog 5.0

Tibase (Locking)



Catalog

Description

IM-LLKEP35

Mist IC L-Link Single 3.5

IM-LLKEP41

Mist IC L-Link Single 4.1

IM-LLKEP50

Mist IC L-Link Single 5.0

Tibase (Non-Locking)



Catalog

Description

IM-3DKEP35

Mist IC 3D Analog 3.5

IM-3DKEP41

Mist IC 3D Analog 4.1

IM-3DKEP50

Mist IC 3D Analog 5.0

Tibase Angulated Screw Channel



Catalog

Description

IM-SLKEP35

Mist IC S-Link 3.5

IM-SLKEP41

Mist IC S-Link 4.1

IM-SLKEP50

Mist IC S-Link 5.0

Ti Blanks 10



Catalog

Description

IM-PF10KEP35

Mist IC PreFit10, 3.5

IM-PF10KEP41

Mist IC PreFit10, 4.1

IM-PF10KEP50

Mist IC PreFit10, 5.0

Ti Blanks 14



Catalog

Description

IM-PF14KEP35

Mist IC PreFit14, 3.5


IM-PF14KEP41

Mist IC PreFit14, 4.1

IM-PF14KEP50

Mist IC PreFit14, 5.0

Final Abutment Screws



Catalog

Description


IM-SCKEP

Mist IC Screw (Pack of 5)

IM-SCKEP UV

Mist IC Angulated Screw (Pack of 5)

Angled Abutment Drivers



Catalog

Description

IM-ALDUVR

Mist IC Angled Latch Driver Reg (13mm)

IM-ALDUVL

Mist IC Angled Latch Driver Long (20mm)

IM-ASDUVR

Mist IC Angled Square Driver Reg (13mm)

IM-ASDUVL

Mist IC Angled Square Driver Long (20mm)

Quad Drivers



Catalog

Description

15272K

Swivel Head, Short

15273K

Swivel Head, Long

45293K

Swivel Head, X-Long (35 mm)

15274K

Latch Type, Short

15275K

Latch Type, Long

45294K

Latch Type, X-Long (35 mm)

15276K

Torque Wrench, Short

15277K

Torque Wrench, Long



Digital Components

Internal TiLobe®



Scan Abutment



Catalog

Description

CORE3D-4805	Scan Abutment Genesis ø3.5
CORE3D-4806	Scan Abutment Genesis ø3.8
CORE3D-4807	Scan Abutment Genesis ø4.5
CORE3D-4808	Scan Abutment Genesis ø5.5
CORE3D-4809	Scan Abutment Genesis ø6.5

Digital Analog



Catalog

Description

CORE3D-3559	Digital Analog PrimaConnex® ø3.5
CORE3D-3560	Digital Analog PrimaConnex® ø4.1
CORE3D-3561	Digital Analog PrimaConnex® ø5.0

Ti Base (Short)



Catalog

Description

45806K	Ti Base Insert ø3.5, Short
45807K	Ti Base Insert ø4.1, Short
45808K	Ti Base Insert ø5.0, Short

Ti Base (Long)



Catalog

Description

45901K	Ti Base Insert ø3.5, Long
45902K	Ti Base Insert ø4.1, Long
45903K	Ti Base Insert ø5.0, Long

Ti Base Pink (Short)



Catalog

Description

45809K	Ti Base Insert ø3.5, Short
45810K	Ti Base Insert ø4.1, Short
45811K	Ti Base Insert ø5.0, Short

Ti Base Pink (Long)



Catalog

Description

45812K	Ti Base Insert ø3.5, Long
45813K	Ti Base Insert ø4.1, Long
45814K	Ti Base Insert ø5.0, Long




In-Office Milling

Internal TiLobe®

Quad Driver
(Torque Wrench)

Quad Driver
(Swivel Head)



Final Abutment Screw



Esthetic Contour
Ti-Base Abutment



Esthetic Contour Ti
Digital Impression Post



Healing Abutment



Ø4.1 Shown

Esthetic Contour Ti Digital Impression Post



Catalog #

45933K
45934K
45935K

Description

Esthetic Contour Ti Digital Impression Post ø3.5
Esthetic Contour Ti Digital Impression Post ø4.1
Esthetic Contour Ti Digital Impression Post ø5.0

Esthetic Contour Ti-Base Abutment



Catalog #

45930K
45931K
45932K

Description


Esthetic Contour Ti-Base Abutment ø3.5
Esthetic Contour Ti-Base Abutment ø4.1
Esthetic Contour Ti-Base Abutment ø5.0

LOCATOR® Abutments



LOCATOR® Abutments

Prosthetic Options



LOCATOR® Abutment Assembly



<i>Catalog #</i>	<i>Description</i>	<i>Cuff Height</i>
45062K	Ø3.5 mm	1.0 mm
45258K	Ø3.5 mm	2.0 mm
45259K	Ø3.5 mm	3.0 mm
45260K	Ø3.5 mm	4.0 mm
45063K	Ø3.5 mm	5.0 mm
45064K	Ø3.5 mm	6.0 mm
45065K	Ø4.1 mm	1.0 mm
45261K	Ø4.1 mm	2.0 mm
45262K	Ø4.1 mm	3.0 mm
45263K	Ø4.1 mm	4.0 mm
45066K	Ø4.1 mm	5.0 mm
45067K	Ø4.1 mm	6.0 mm
45068K	Ø5.0 mm	1.0 mm
45264K	Ø5.0 mm	2.0 mm
45265K	Ø5.0 mm	3.0 mm
45069K	Ø5.0 mm	4.0 mm
45070K	Ø5.0 mm	5.0 mm
45071K	Ø5.0 mm	6.0 mm

LOCATOR® Abutment Accessories



<i>Catalog #</i>	<i>Description</i>
L8519-2	Denture Cap with Black Processing Male, 2-Pack Clear, Pink, Blue Retention, Block out spacers
45266K	LOCATOR® Replacement Male, Clear, Standard Retention, 4-Pack
45267K	LOCATOR® Replacement Male, Pink, Light Retention, 4-Pack
45268K	LOCATOR® Replacement Male, Blue, (ELR) Extra Light Retention, 4-Pack
45269K	LOCATOR® Replacement Male, Green, (ER) Extended Range, 4-Pack
45270K	LOCATOR® Female Analog, 4-Pack
45271K	LOCATOR® Impression Coping, 4-Pack
45273K	LOCATOR® Core Tool
45274K	LOCATOR® Torque Wrench Driver



LOCATOR® R-Tx™ Abutment

Prosthetic Options

LOCATOR R-Tx™ Abutment



Catalog

Description

Cuff Height

LOCR-TX 31200-01
LOCR-TX 31200-02
LOCR-TX 31200-03
LOCR-TX 31200-04
LOCR-TX 31200-05
LOCR-TX 31200-06

Ø3.5 mm
Ø3.5 mm
Ø3.5 mm
Ø3.5 mm
Ø3.5 mm
Ø3.5 mm

1.0 mm
2.0 mm
3.0 mm
4.0 mm
5.0 mm
6.0 mm

LOCR-TX 31201-01
LOCR-TX 31201-02
LOCR-TX 31201-03
LOCR-TX 31201-04
LOCR-TX 31201-05
LOCR-TX 31201-06

Ø4.1 mm
Ø4.1 mm
Ø4.1 mm
Ø4.1 mm
Ø4.1 mm
Ø4.1 mm

1.0 mm
2.0 mm
3.0 mm
4.0 mm
5.0 mm
6.0 mm

LOCR-TX 31202-01
LOCR-TX 31202-02
LOCR-TX 31202-03
LOCR-TX 31202-04
LOCR-TX 31202-05
LOCR-TX 31202-06

Ø5.0 mm
Ø5.0 mm
Ø5.0 mm
Ø5.0 mm
Ø5.0 mm
Ø5.0 mm

1.0 mm
2.0 mm
3.0 mm
4.0 mm
5.0 mm
6.0 mm

LOCATOR R-Tx™ Abutment Accessories



Catalog

Description

LOCR-TX 30013-01

Denture Attachment Processing Assembly, 4-Pack

LOCR-TX 30001-01

Zero Retention Insert, 4-Pack

LOCR-TX 30002-01

Low Retention Insert, 4-Pack

LOCR-TX 30003-01

Medium Retention Insert, 4-Pack

LOCR-TX 30004-01

High Retention Insert, 4-Pack

LOCR-TX 30014-01

3.35 mm Abutment Analog, 4-Pack

LOCR-TX 30015-01

4.0 mm Abutment Analog, 4-Pack

LOCR-TX 30016-01

5.0 mm Abutment Analog, 4-Pack

LOCR-TX 30017-01

Impression Coping, 4-Pack

LOCR-TX 30021-01

Retention Insert Tool

LOCR-TX 04936

.050"/1.25 mm Hex, Torque Wrench

LOCR-TX 08007

.050"/1.25 mm Hex Driver with Thumb Knob

LOCR-TX 08008

.050"/1.25 mm Hex, Latch Lock Driver

LOCR-TX 30012-01

Processing Insert, 4-Pack

LOCR-TX 30018-01

Processing Spacer, 4-Pack

LOCATOR F-Tx®
Abutment



Catalog #

Description

Cuff Height

LOCF-TX 11200-01	Ø3.5 mm	1.0 mm
LOCF-TX 11200-02	Ø3.5 mm	2.0 mm
LOCF-TX 11200-03	Ø3.5 mm	3.0 mm
LOCF-TX 11200-04	Ø3.5 mm	4.0 mm
LOCF-TX 11200-05	Ø3.5 mm	5.0 mm
LOCF-TX 11201-01	Ø4.1 mm	1.0 mm
LOCF-TX 11201-02	Ø4.1 mm	2.0 mm
LOCF-TX 11201-03	Ø4.1 mm	3.0 mm
LOCF-TX 11201-04	Ø4.1 mm	4.0 mm
LOCF-TX 11201-05	Ø4.1 mm	5.0 mm
LOCF-TX 11202-01	Ø5.0 mm	1.0 mm
LOCF-TX 11202-02	Ø5.0 mm	2.0 mm
LOCF-TX 11202-03	Ø5.0 mm	3.0 mm
LOCF-TX 11202-04	Ø5.0 mm	4.0 mm
LOCF-TX 11202-05	Ø5.0 mm	5.0 mm

LOCATOR F-Tx™
Abutment Accessories




Catalog #

Description

LOCF-TX 10028-01	Manual Driver
LOCF-TX 10029-01	Latch Driver
LOCF-TX 10030-01	Processing Cap, 2Pk
LOCF-TX 10032-01	Processing Ball, 2Pk
LOCF-TX 10033-01	Low Retention Ball, 2Pk
LOCF-TX 10034-01	Medium Retention Ball, 2Pk
LOCF-TX 10035-01	High Retention Ball, 2 Pk
LOCF-TX 10036-01	Healing Cap, 2Pk
LOCF-TX 10037-01	Abutment Analog, 2Pk
LOCF-TX 10039-01	Polishing Cap, 2Pk
LOCF-TX 10040-01	Waxing Cap, 2Pk
LOCF-TX 10041-01	Denture Attachment Housing, 2Pk
LOCF-TX 10042-01	Retention Ball Hex Driver
LOCF-TX 10043-01	Impression Coping, 2Pk
LOCF-TX 10094	Tip Wrench
LOCF-TX 10095	Seating Tip
LOCF-TX 10096	Removal Tip with Threaded Adapter
LOCF-TX 10097	Replacement Wire
LOCF-TX 10098	Threaded Adapter
LOCF-TX 10099	Wire and Lever Wrench
LOCF-TX 100100	Seating and Removal Tool

Surgical Kit & Instruments



Catalog #	Surgical Options
G31000	Surgical Cassette (empty)
G31001	A. Genesis Surgical Kit for Ø3.5-Ø4.5 mm
G31002	B. Genesis Drill Pack for Ø5.5/Ø6.5 mm

Catalog #	Description	Catalog #	Description
G31001	A. Genesis Surgical Kit for Ø3.5-Ø4.5 mm <i>Includes:</i>	G31036	Implant Depth Gauge 4.5 x 8.5-10 mm
G31000	Surgical Cassette	G31037	Implant Depth Gauge 4.5 x 11.5-13 mm
G21244	2.0 mm Full Depth Twist Drill	G31038	Implant Depth Gauge 4.5 x 16 mm
G21211	Parallel Pin 4-Pack	G21290	3.5/3.8 Implant Driver, Latch, Short
G31008	Final Tapered Drill 3.5 x 8.5 mm	G21291	3.5/3.8 Implant Driver, Latch, Long
G31009	Final Tapered Drill 3.5 x 10 mm	G21292	4.5 Implant Driver, Latch, Short
G31010	Final Tapered Drill 3.5 x 11.5 mm	G21293	4.5 Implant Driver, Latch, Long
G31011	Final Tapered Drill 3.5 x 13 mm	G21274	3.5/3.8 Implant Driver, Ratchet
G31012	Final Tapered Drill 3.5 x 16 mm	G21275	4.5 Implant Driver, Ratchet
G31013	Final Tapered Drill 3.8 x 8.5 mm	K15695	Surgical Ratchet/Torque Wrench
G31014	Final Tapered Drill 3.8 x 10 mm	15224K	Ratchet Adapter
G31015	Final Tapered Drill 3.8 x 11.5 mm	G71612	Quad Driver, Torque Wrench, Short
G31016	Final Tapered Drill 3.8 x 13 mm	G71613	Quad Driver, Torque Wrench, Long
G31017	Final Tapered Drill 3.8 x 16 mm	G71608	Quad Driver, Swivel Head, Short
G31018	Final Tapered Drill 4.5 x 8.5 mm	G71609	Quad Driver, Swivel Head, Long
G31019	Final Tapered Drill 4.5 x 10 mm	G21204	Drill Extender
G31020	Final Tapered Drill 4.5 x 11.5 mm	G21277	Conversion Handle, Latch Type
G31021	Final Tapered Drill 4.5 x 13 mm	G21278	Conversion Handle, Square Type
G31022	Final Tapered Drill 4.5 x 16 mm	100001-EN01	Surgical Manual
G31031	Implant Depth Gauge 3.5 x 10-11.5 mm	100009-EN01	Drill Usage Chart
G31032	Implant Depth Gauge 3.5 x 13-16 mm	100005-EN01	Wall Chart
G31033	Implant Depth Gauge 3.8 x 8.5-10 mm	100027-EN01	Implant Driver, Latch IFU
G31034	Implant Depth Gauge 3.8 x 11.5 -13 mm	MK40091	Genesis Folder
G31035	Implant Depth Gauge 3.8 x 16 mm		



<i>Catalog #</i>	<i>Description</i>
G31002	B. Genesis Drill Pack for Ø5.5/6.5 mm <i>Includes:</i>
G31023	Final Tapered Drill 5.5 x 8.5 mm
G31024	Final Tapered Drill 5.5 x 10 mm
G31025	Final Tapered Drill 5.5 x 11.5 mm
G31026	Final Tapered Drill 5.5 x 13 mm
G31027	Final Tapered Drill 6.5 x 8.5 mm
G31028	Final Tapered Drill 6.5 x 10 mm
G31029	Final Tapered Drill 6.5 x 11.5 mm
G31030	Final Tapered Drill 6.5 x 13 mm
G31039	Implant Depth Gauge 5.5 x 8.5-10 mm
G31040	Implant Depth Gauge 5.5 x 11.5-13 mm
G31041	Implant Depth Gauge 6.5 x 8.5-10 mm
G31042	Implant Depth Gauge 6.5 x 11.5-13 mm
G21294	5.5/6.5 Implant Driver, Latch, Short
G21295	5.5/6.5 Implant Driver, Latch, Long
G21276	5.5/6.5 Implant Driver, Ratchet
100027-EN01	Implant Driver, Latch IFU

Initial Drill



<i>Catalog #</i>	<i>Description</i>
G21243	Spade Drill
G21244	2.0 mm Full Depth Twist Drill


Final Tapered Drills



<i>Catalog #</i>	<i>Description</i>
G31008	Ø3.5 x 8.5 mm
G31009	Ø3.5 x 10 mm
G31010	Ø3.5 x 11.5 mm
G31011	Ø3.5 x 13 mm
G31012	Ø3.5 x 16 mm
G31013	Ø3.8 x 8.5 mm
G31014	Ø3.8 x 10 mm
G31015	Ø3.8 x 11.5 mm
G31016	Ø3.8 x 13 mm
G31017	Ø3.8 x 16 mm
G31018	Ø4.5 x 8.5 mm
G31019	Ø4.5 x 10 mm
G31020	Ø4.5 x 11.5 mm
G31021	Ø4.5 x 13 mm
G31022	Ø4.5 x 16 mm
G31023	Ø5.5 x 8.5 mm
G31024	Ø5.5 x 10 mm
G31025	Ø5.5 x 11.5 mm
G31026	Ø5.5 x 13 mm



Final Tapered Drills (cont.)



Catalog

G31027
G31028
G31029
G31030

Description

Ø6.5 x 8.5 mm
Ø6.5 x 10 mm
Ø6.5 x 11.5 mm
Ø6.5 x 13 mm

Depth Gauges



Catalog

G31031
G31032

Description

Ø3.5 x 10-11.5 mm
Ø3.5 x 13-16 mm

G31033
G31034
G31035

Ø3.8 x 8.5-10 mm
Ø3.8 x 11.5-13 mm
Ø3.8 x 16 mm

G31036
G31037
G31038

Ø4.5 x 8.5-10 mm
Ø4.5 x 11.5-13 mm
Ø4.5 x 16 mm

G31039
G31040

Ø5.5 x 8.5-10 mm
Ø5.5 x 11.5-13 mm

G31041
G31042

Ø6.5 x 8.5-10 mm
Ø6.5 x 11.5-13 mm

Surgical Taps



Catalog

G31045
G31046

Description

Ø3.5 x 10-11.5 mm
Ø3.5 x 13-15 mm

G21008

Ø3.8 x 11.5-16 mm

G21009
G21010

Ø4.5 x 10-13 mm
Ø4.5 x 16 mm

G21011
G21012

Ø5.5 x 8.5 mm
Ø5.5 x 10-13 mm

G21013
G21014

Ø6.5 x 8.5 mm
Ø6.5 x 10-13 mm

Parallel Pin



Catalog

G21211

Description

4-Pack

Drill Extender



Catalog

G21204



Implant Driver, Latch



Catalog

G21290
G21291

G21292
G21293

G21294
G21295

Description

Ø3.5/3.8 mm, Short
Ø3.5/3.8 mm, Long

Ø4.5 mm, Short
Ø4.5 mm, Long

Ø5.5/6.5 mm, Short
Ø5.5/6.5 mm, Long

Implant Driver, Ratchet



Catalog

G21274
G21275
15284K

Description

Ø3.5/3.8 mm
Ø4.5 mm
Ø5.5/6.5 mm

Surgical Ratchet/ Torque Wrench



Catalog

K15695

Description

Surgical Ratchet/Torque Wrench

Ratchet Adapter



Catalog

15224K

Conversion Handles




Catalog

G21277
G21278

Description

Conversion Handle, Latch Type
Conversion Handle, Square Type

Quad Drivers



Catalog

G71608
G71609
45293K

G71610
G71611
45294K

G71612
G71613

Description

Swivel Head, Short
Swivel Head, Long
Swivel Head, X-Long (35 mm)

Latch Type, Short
Latch Type, Long
Latch Type, X-Long (35 mm)

Torque Wrench, Short
Torque Wrench, Long

BESTELLSCHEIN

FAX: 02064 625 95 80 oder 0800 078 63 27

Sie erreichen unsere Bestellannahme
telefonisch Mo. - Do. 7:45 - 18:00 Uhr
Fr. 7:45 - 15:00 Uhr
Fax / E-Mail rund um die Uhr

RUNDAS GmbH
Amalienstr. 62
46537 Dinslaken
Tel. +49 2064 625 95 50
Fax +49 2064 625 95 80 oder 0800 078 63 27 (kostenfrei)

Absender

Name:

Kunden-Nr.:

Telefon/Fax:

E-Mail:

Datum/

Unterschrift:

☐ Bitte senden Sie mir Ihren aktuellen Katalog.

Artikelnummer	Artikelbezeichnung	Menge	Einzelpreis Euro *

* zzgl. MwSt. und Versand

Ab 100,- € Netto-Warenwert Lieferung frei Haus!
Unter 100,- € Netto-Warenwert berechnen wir pauschal 5,- € Versandkosten.
Alle Preise sind netto zzgl. gesetzl. MwSt. und verstehen sich als Stückpreise.
Es gelten unsere AGB. Zahlungsziel: 10 Tage – 2 % Skonto, 30 Tage – netto Kasse

RUNDAS GmbH • Amalienstr. 62 • 46537 Dinslaken
Tel. +49 2064 625 95 50
Fax +49 2064 625 95 80 oder 0800 078 63 27 (kostenfrei)
info@rundas.de • www.rundas.de



**Keystone Dental, Inc.
Global Headquarters**

154 Middlesex Turnpike
Burlington, MA 01803 USA
General Inquiries: 781-328-3600
Customer Relations: 866-902-9272 or 781-328-3490
Fax: 866-903-9272 or 781-328-3400
Email: info@keystonedental.com

Keystone Dental, Inc.

5 Holland, Bldg 209
Irvine, CA 92618 USA
Customer Service Number: 866-902-9272
Fax: 866-903-9272 or 781-328-3400
Email: info@keystonedental.com



EMERGO Europe B.V.



Prinsessegracht 20
2514 AP The Hague
The Netherlands

Vertrieb / Distribution

Deutschland, Belgien,
Niederlande, Luxemburg,
Frankreich:

RUNDAS GmbH
Amalienstr. 62
46537 Dinslaken

Tel.: 02064 625 95 50
Fax: 02064 625 95 80

E-Mail: info@rundas.de
Internet: www.rundas.de



Caution, consult accompanying documents

Trademark Acknowledgements

TiLobe®, PrimaPlus, Quick-Cap and "Smarter Thinking. Simpler Design." are trademarks and registered trademarks of Keystone Dental, Inc.