

Smarter Thinking. Simpler Design.



PrimaConnex[®]

Product Catalog



PrimaConnex®

Ø3.5 mm PrimaConnex® Implants	6
Ø3.5 mm Prosthetic Options	7
Ø3.5 mm Prosthetic Components	8
Ø4.1 mm PrimaConnex® Implants	12
Ø4.1 mm Prosthetic Options	13
Ø4.1 mm Prosthetic Components	14
Ø5.0 mm PrimaConnex® Implants	18
Ø5.0 mm Prosthetic Options	19
Ø5.0 mm Prosthetic Components	20
LOCATOR® Prosthetic Options	25
LOCATOR R-Tx™ Abutment	26
LOCATOR® Abutment	27
PrimaConnex® Surgical Kit	29
PrimaConnex® Guided Surgical Kit.....	35
Terms & Conditions	38


Product specifications are subject to change without notice.
 Items illustrated are not to scale.
 Certain products may not be available in all countries.



PrimaConnex[®]



PrimaConnex® Implants Internal TiLobe




PrimaConnex® Implants Tapered*

Catalog

Implant Size

15413K	Ø3.5 x 10 mm Length
15414K	Ø3.5 x 11.5 mm Length
15415K	Ø3.5 x 13 mm Length
15416K	Ø3.5 x 15 mm Length
15613K	Ø3.5 x 10 mm Length TC
15614K	Ø3.5 x 11.5 mm Length TC
15615K	Ø3.5 x 13 mm Length TC
15616K	Ø3.5 x 15 mm Length TC



PrimaConnex® Implants Straight*

Catalog

Implant Size

15425K	Ø3.3 x 10 mm Length
15426K	Ø3.3 x 11.5 mm Length
15427K	Ø3.3 x 13 mm Length
15428K	Ø3.3 x 15 mm Length
15625K	Ø3.3 x 10 mm Length TC
15626K	Ø3.3 x 11.5 mm Length TC
15627K	Ø3.3 x 13 mm Length TC
15628K	Ø3.3 x 15 mm Length TC



Cover Screw

Catalog

Description

15300K	Ø3.5 mm
---------------	---------



Healing Abutments

Catalog

Description


Flare

Cuff Height

15301K	Ø3.5 mm	Straight	3.0 mm
15302K	Ø3.5 mm	Straight	5.0 mm
15303K	Ø3.5 mm	Straight	7.0 mm
15317K	Ø3.5 mm	4.0 mm	3.0 mm
15318K	Ø3.5 mm	4.0 mm	5.0 mm
15304K	Ø3.5 mm	5.0 mm	3.0 mm
15305K	Ø3.5 mm	5.0 mm	5.0 mm

*All PrimaConnex® implant sizes are available with either a 1.0 mm machined collar or Total Coverage (TC).

The Ø3.5 mm PrimaConnex® can be restored with straight, Ø4.0 and Ø5.0 mm flared prosthetic components. An Abutment Screw is included with all TiLobe abutments.





Prosthetic Components

Internal TiLobe

Open & Closed Tray Impression Post



Catalog

Description

Flare

45130K
45136K
45131K

Ø3.5 mm
Ø3.5 mm
Ø3.5 mm

Straight
4.0 mm
5.0 mm

Open Tray Impression Post



Catalog

Description

Flare

45904K
45906K

Ø3.5/3.8 mm
Ø3.5/3.8 mm

4.0 mm
5.0 mm

Closed Tray Impression Post



Catalog

Description

Flare

45905K
45907K

Ø3.5/3.8 mm
Ø3.5/3.8 mm

4.0 mm
5.0 mm

Implant Analog



Catalog

Description

45140K
45143K

Ø3.5 mm, 1-Pack
Ø3.5 mm, 10-Pack

Temporary Abutment



Catalog

Description

Flare

45118K
45119K

45112K
45113K

Ø3.5 mm, Locking
Ø3.5 mm, Non-Locking

Ø3.5 mm, Locking
Ø3.5 mm, Non-Locking

4.0 mm
4.0 mm

5.0 mm
5.0 mm

Ti Temporary Abutment



Catalog

Description

45075K
45076K

Ø3.5 mm, Locking
Ø3.5 mm, Non-Locking

Esthetic Contour Ti Abutment



Cuff Height
(Buccal)

Catalog

Description

Flare

Cuff Height

45046K
45047K
45048K
45049K

Ø3.5 mm, Straight
Ø3.5 mm, Straight
Ø3.5 mm, 15° Angled
Ø3.5 mm, 15° Angled

4.0 mm
4.0 mm
4.0 mm
4.0 mm

1.0 mm
2.0 mm
1.0 mm
2.0 mm

45030K
45031K
45032K
45033K

Ø3.5 mm, Straight
Ø3.5 mm, Straight
Ø3.5 mm, 15° Angled
Ø3.5 mm, 15° Angled

5.0 mm
5.0 mm
5.0 mm
5.0 mm

1.0 mm
3.0 mm
1.0 mm
3.0 mm

UCLA Abutment



Catalog

Description

Cuff Height

45100K
45101K

Gold/Plastic Sleeve, Locking
Gold/Plastic Sleeve, Non-Locking

1.0 mm
1.0 mm



**Replacement
Final Abutment Screw**

Catalog #

Description

45060K

1-Pack



Laboratory Screw

Catalog #

Description

45061K

10-Pack



Quick-Packs

Catalog #

Description

Flare

Cuff Height

45400K

Ø3.5 mm

5.0 mm

1.0 mm

45401K

Ø3.5 mm

5.0 mm

2.0 mm

45402K

Ø3.5 mm

5.0 mm


3.0 mm

45403K

Ø3.5 mm

5.0 mm

4.0 mm



Quick-Abutment

Catalog #

Description

Flare

Cuff Height

45420K

Ø3.5 mm

5.0 mm

1.0 mm

45421K

Ø3.5 mm

5.0 mm

2.0 mm

45422K

Ø3.5 mm

5.0 mm

3.0 mm

45423K

Ø3.5 mm

5.0 mm

4.0 mm



Quick-Abutment Quick-Cap®

Catalog #

Description

Flare

45432K

Ø3.5 mm, 1-Pack

5.0 mm



**Quick-Abutment
Temporary/Healing Cap**

Catalog #

Description

Flare

45209K

Ø3.5 mm

5.0 mm



Quick-Abutment Analog

Catalog #

Description

Flare

45445K

Ø3.5 mm

5.0 mm



**Quick-Abutment
Plastic Sleeve**

Catalog #

Description

Flare

45215K

Ø3.5 mm, Locking

5.0 mm

45218K

Ø3.5 mm, Non-Locking

5.0 mm



**Quick-Abutment
Ceramic Coping**

Catalog #

Description

Flare

45301K

Ø3.5 mm

5.0 mm



Reamer Tool for Quick-Abutment Plastic Sleeve

Catalog #

Description

45311K

Reamer Tool and Centering Pins




Prepable Abutment

Catalog #

Description

45042K

Ø3.5 mm




Ti-Base Abutment (Short)

Catalog #

Description

45806K

Ø3.5 mm




Pink Ti-Base Abutment (Short)

Catalog #

Description

45007K

Ø3.5/3.8 mm



Ti-Base Abutment (Long)

Catalog #

Description

45901K


Ø3.5 mm

Multi Unit Abutment

Catalog #

Description

Cuff Height



45312K

Ø3.5/3.8 mm, Straight

1.0 mm

45313K

Ø3.5/3.8 mm, Straight

2.0 mm

45314K

Ø3.5/3.8 mm, Straight

3.0 mm

45321K

Ø3.5/3.8 mm, Angled, 17°

3.0 mm

45322K

Ø3.5/3.8 mm, Angled, 17°

4.0 mm

45936K

Ø3.5/3.8 mm, Angled, 30°

3.0 mm

45937K

Ø3.5/3.8 mm, Angled, 30°

4.0 mm

45938K

Ø3.5/3.8 mm, Angled, 30°

5.0 mm

**Multi Unit Abutment
Healing Cap**



Catalog #

Description

45330K

Ø3.5/3.8 mm

**Multi Unit Abutment
Impression Post**



Catalog #

Description

**45327K
45328K**

Open Tray, Ø3.5/3.8 mm with screw
Closed Tray, Ø3.5/3.8 mm

**Multi Unit Abutment
Analog**



Catalog #

Description

45329K

Ø3.5/3.8 mm

**Multi-Unit Abutment
Ti Temporary Sleeve**



Catalog #

Description

45331K

Ø3.5/3.8 mm

**Multi Unit Abutment
Plastic Sleeve**



Catalog #

Description

45332K

Ø3.5/3.8 mm

**Multi Unit Abutment
Prosthetic Screw**



Catalog #

Description

45333K

1-Pack

**Replacement Abutment
Screw for Multi Unit
Angled Abutment**




Catalog #

Description

45334K

1-Pack




**PrimaConnex® Implants
Tapered***

Catalog #

Implant Size

15417K	Ø4.1 x 10 mm Length
15418K	Ø4.1 x 11.5 mm Length
15419K	Ø4.1 x 13 mm Length
15420K	Ø4.1 x 15 mm Length
15617K	Ø4.1 x 10 mm Length TC
15618K	Ø4.1 x 11.5 mm Length TC
15619K	Ø4.1 x 13 mm Length TC
15620K	Ø4.1 x 15 mm Length TC



**PrimaConnex® Implants
Straight***

Catalog #

Implant Size

15429K	Ø4.0 x 8 mm Length
15430K	Ø4.0 x 10 mm Length
15431K	Ø4.0 x 11.5 mm Length
15432K	Ø4.0 x 13 mm Length
15433K	Ø4.0 x 15 mm Length
15629K	Ø4.0 x 8 mm Length TC
15630K	Ø4.0 x 10 mm Length TC
15631K	Ø4.0 x 11.5 mm Length TC
15632K	Ø4.0 x 13 mm Length TC
15633K	Ø4.0 x 15 mm Length TC



Cover Screw

Catalog #

Description

15306K	Ø4.1 mm
---------------	---------



Healing Abutments

Catalog #

Description


Flare

Cuff Height

15307K	Ø4.1 mm	Straight	3.0 mm
15308K	Ø4.1 mm	Straight	5.0 mm
15309K	Ø4.1 mm	Straight	7.0 mm
15319K	Ø4.1 mm	5.0 mm	3.0 mm
15320K	Ø4.1 mm	5.0 mm	5.0 mm
15310K	Ø4.1 mm	6.0 mm	3.0 mm
15311K	Ø4.1 mm	6.0 mm	5.0 mm

*All PrimaConnex® implant sizes are available with either a 1.0 mm machined collar or Total Coverage (TC).

The Ø4.1 mm PrimaConnex® can be restored with straight, Ø5.0 and Ø6.0 mm flared prosthetic components. An Abutment Screw is included with all TiLobe abutments.



**Open & Closed Tray
Impression Post***Catalog #**Description**Flare***45132K**

Ø4.1 mm

Straight

45137K

Ø4.1 mm

5.0 mm

45133K

Ø4.1 mm

6.0 mm

**Open Tray
Impression Post***Catalog #**Description**Flare***45908K**

Ø4.1/4.5 mm

5.0 mm

45910K

Ø4.1/4.5 mm

6.0 mm

**Closed Tray
Impression Post***Catalog #**Description**Flare***45909K**

Ø4.1/4.5 mm

5.0 mm

45911K

Ø4.1/4.5 mm

6.0 mm

Implant Analog*Catalog #**Description***45141K**

Ø4.1 mm, 1-Pack

45144K

Ø4.1 mm, 10-Pack

Temporary Abutment*Catalog #**Description**Flare***45120K**

Ø4.1 mm, Locking

5.0 mm

45121K

Ø4.1 mm, Non-Locking

5.0 mm

45114K

Ø4.1 mm, Locking

6.0 mm

45115K

Ø4.1 mm, Non-Locking

6.0 mm

Ti Temporary Abutment*Catalog #**Description***45077K**

Ø4.1 mm, Locking

45078K

Ø4.1 mm, Non-Locking

**Esthetic Contour Ti
Abutment**Cuff Height —
(Buccal)*Catalog #**Description**Flare**Cuff Height***45050K**

Ø4.1 mm, Straight

5.0 mm

0.5 mm

45051K

Ø4.1 mm, Straight

5.0 mm

1.5 mm

45052K

Ø4.1 mm, 15° Angled

5.0 mm

0.5 mm

45053K

Ø4.1 mm, 15° Angled

5.0 mm

1.5 mm

45034K

Ø4.1 mm, Straight

6.0 mm

1.0 mm

45035K

Ø4.1 mm, Straight

6.0 mm

3.0 mm

45036K

Ø4.1 mm, 15° Angled

6.0 mm

1.0 mm

45037K

Ø4.1 mm, 15° Angled

6.0 mm

3.0 mm

UCLA Abutment*Catalog #**Description**Cuff Height***45102K**

Gold/Plastic Sleeve, Locking

1.0 mm

45103K

Gold/Plastic Sleeve, Non-Locking

1.0 mm



**Replacement
Final Abutment Screw**

<i>Catalog #</i>	<i>Description</i>
45060K	1-Pack




Laboratory Screw

<i>Catalog #</i>	<i>Description</i>
45061K	10-Pack



Quick-Packs

<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
45412K	Ø4.1 mm	5.0 mm	0.5 mm
45413K	Ø4.1 mm	5.0 mm	1.5 mm
45414K	Ø4.1 mm	5.0 mm	2.5 mm
45415K	Ø4.1 mm	5.0 mm	4.0 mm
45404K	Ø4.1 mm	6.0 mm	0.5 mm
45405K	Ø4.1 mm	6.0 mm	1.5 mm
45406K	Ø4.1 mm	6.0 mm	2.5 mm
45407K	Ø4.1 mm	6.0 mm	4.0 mm



Quick-Abutment

<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
45435K	Ø4.1 mm	5.0 mm	0.5 mm
45436K	Ø4.1 mm	5.0 mm	1.5 mm
45437K	Ø4.1 mm	5.0 mm	2.5 mm
45438K	Ø4.1 mm	5.0 mm	4.0 mm
45424K	Ø4.1 mm	6.0 mm	0.5 mm
45425K	Ø4.1 mm	6.0 mm	1.5 mm
45426K	Ø4.1 mm	6.0 mm	2.5 mm
45427K	Ø4.1 mm	6.0 mm	4.0 mm



**Quick-Abutment
Quick-Cap®**

<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
45443K	Ø4.1 mm, 1-Pack	5.0 mm
45433K	Ø4.1 mm, 1-Pack	6.0 mm



**Quick-Abutment
Temporary/Healing Cap**

<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
45236K	Ø4.1 mm	5.0 mm
45210K	Ø4.1 mm	6.0 mm



**Quick-Abutment
Analog**

<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
45448K	Ø4.1 mm	5.0 mm
45446K	Ø4.1 mm	6.0 mm

**Quick-Abutment
Plastic Sleeve***Catalog #**Description**Flare***45240K**

Ø4.1 mm, Locking

5.0 mm

45242K

Ø4.1 mm, Non-Locking

5.0 mm

45216K

Ø4.1 mm, Locking

6.0 mm

45219K

Ø4.1 mm, Non-Locking

6.0 mm

**Quick-Abutment
Ceramic Coping***Catalog #**Description**Flare***45244K**

Ø4.1 mm

5.0 mm

45302K

Ø4.1 mm

6.0 mm

**Reamer Tool for Quick-
Abutment Plastic Sleeve***Catalog #**Description***45311K**


Reamer Tool and Centering Pins

**Prepable Abutment***Catalog #**Description***45043K**

Ø4.1 mm

**Ti-Base Abutment (Short)***Catalog #**Description***45807K**

Ø4.1 mm


**Pink Ti-Base
Abutment (Short)***Catalog #**Description***45008K**

Ø4.1/4.5 mm

**Ti-Base Abutment (Long)***Catalog #**Description***45902K**

Ø4.1 mm

Multi Unit Abutment



Catalog

Description

Cuff Height

45315K	Ø4.1/4.5 mm, Straight	1.0 mm
45316K	Ø4.1/4.5 mm, Straight	2.0 mm
45317K	Ø4.1/4.5 mm, Straight	3.0 mm
45318K	Ø4.1/4.5 mm, Straight	4.0 mm
45323K	Ø4.1/4.5 mm, Angled, 17°	3.0 mm
45324K	Ø4.1/4.5 mm, Angled, 17°	4.0 mm
45325K	Ø4.1/4.5 mm, Angled, 30°	4.0 mm
45326K	Ø4.1/4.5 mm, Angled, 30°	5.0 mm

Multi Unit Abutment Healing Cap



Catalog

Description

45330K	Ø4.1/4.5 mm
---------------	-------------

Multi Unit Abutment Impression Post



Catalog

Description

45327K	Open Tray, Ø4.1/4.5 mm with screw
45328K	Closed Tray, Ø4.1/4.5 mm

Multi Unit Abutment Analog



Catalog

Description

45329K	Ø4.1/4.5 mm
---------------	-------------

Multi-Unit Abutment Ti Temporary Sleeve



Catalog

Description

45331K	Ø4.1/4.5 mm
---------------	-------------

Multi Unit Abutment Plastic Sleeve



Catalog

Description

45332K	Ø4.1/4.5 mm
---------------	-------------

Multi Unit Abutment Prosthetic Screw



Catalog

Description

45333K	1-Pack
---------------	--------

Replacement Abutment Screw for Multi Unit Angled Abutment




Catalog

Description

45334K	1-Pack
---------------	--------



PrimaConnex® Implants Internal TiLobe




PrimaConnex® Implants Tapered*

Catalog

Implant Size

15421K	Ø5.0 x 10 mm Length
15422K	Ø5.0 x 11.5 mm Length
15423K	Ø5.0 x 13 mm Length
15424K	Ø5.0 x 15 mm Length
15621K	Ø5.0 x 10 mm Length TC
15622K	Ø5.0 x 11.5 mm Length TC
15623K	Ø5.0 x 13 mm Length TC
15624K	Ø5.0 x 15 mm Length TC



PrimaConnex® Implants Straight*

Catalog

Implant Size

15434K	Ø5.0 x 8 mm Length
15435K	Ø5.0 x 10 mm Length
15436K	Ø5.0 x 11.5 mm Length
15437K	Ø5.0 x 13 mm Length
15634K	Ø5.0 x 8 mm Length TC
15635K	Ø5.0 x 10 mm Length TC
15636K	Ø5.0 x 11.5 mm Length TC
15637K	Ø5.0 x 13 mm Length TC



Cover Screw

Catalog

Description

15312K	Ø5.0 mm
---------------	---------



Healing Abutments

Catalog

Description


Flare

Cuff Height

15313K	Ø5.0 mm	Straight	3.0 mm
15314K	Ø5.0 mm	Straight	5.0 mm
15321K	Ø5.0 mm	5.0 mm	3.0 mm
15322K	Ø5.0 mm	5.0 mm	5.0 mm
15315K	Ø5.0 mm	6.0 mm	3.0 mm
15316K	Ø5.0 mm	6.0 mm	5.0 mm

*All PrimaConnex® implant sizes are available with either a 1.0 mm machined collar or Total Coverage (TC).

The Ø5.0 mm PrimaConnex® can be restored with straight, Ø5.0 and Ø6.0 mm flared prosthetic components. An Abutment Screw is included with all TiLobe abutments.



**Open & Closed Tray
Impression Post***Catalog #**Description**Flare***45134K**

Ø5.0 mm

Straight

45138K

Ø5.0 mm

5.0 mm

45135K

Ø5.0 mm

6.0 mm

**Open Tray
Impression Post***Catalog #**Description**Flare***45916K**

Ø5.5/6.5 mm

7.0 mm

45918K

Ø5.5/6.5 mm

8.0 mm

**Closed Tray
Impression Post***Catalog #**Description**Flare***45917K**

Ø5.5/6.5 mm

7.0 mm

45919K

Ø5.5/6.5 mm

8.0 mm

**Implant Analog***Catalog #**Description***45142K**

Ø5.0 mm, 1-Pack

45145K

Ø5.0 mm, 10-Pack

**Temporary Abutment***Catalog #**Description**Flare***45122K**

Ø5.0 mm, Locking

5.0 mm

45123K

Ø5.0 mm, Non-Locking

5.0 mm

45116K

Ø5.0 mm, Locking

6.0 mm

45117K

Ø5.0 mm, Non-Locking

6.0 mm

**Ti Temporary Abutment***Catalog #**Description***45079K**

Ø5.0 mm, Locking

45080K

Ø5.0 mm, Non-Locking

Cuff Height
(Buccal)**Esthetic Contour Ti
Abutment***Catalog #**Description**Flare**Cuff Height***45054K**

Ø5.0 mm, Straight

5.0 mm

0.5 mm

45055K

Ø5.0 mm, Straight

5.0 mm

1.5 mm

45056K

Ø5.0 mm, 15° Angled

5.0 mm

0.5 mm

45057K

Ø5.0 mm, 15° Angled

5.0 mm

1.5 mm

45038K

Ø5.0 mm, Straight

6.0 mm

1.0 mm

45039K

Ø5.0 mm, Straight

6.0 mm

3.0 mm

45040K

Ø5.0 mm, 15° Angled

6.0 mm

1.0 mm

**UCLA Abutment***Catalog #**Description**Cuff Height***45104K**

Gold/Plastic Sleeve, Locking

1.0 mm

45105K

Gold/Plastic Sleeve, Non-Locking

1.0 mm



**Replacement
Final Abutment Screw**

<i>Catalog #</i>	<i>Description</i>
45060K	1-Pack



Laboratory Screw


<i>Catalog #</i>	<i>Description</i>
45061K	10-Pack

Quick-Packs



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
45416K	Ø5.0 mm	5.0 mm	0.5 mm
45417K	Ø5.0 mm	5.0 mm	1.5 mm
45418K	Ø5.0 mm	5.0 mm	2.5 mm
45419K	Ø5.0 mm	5.0 mm	4.0 mm
45408K	Ø5.0 mm	6.0 mm	0.5 mm
45409K	Ø5.0 mm	6.0 mm	1.5 mm
45410K	Ø5.0 mm	6.0 mm	2.5 mm
45411K	Ø5.0 mm	6.0 mm	4.0 mm

Quick-Abutment



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
45439K	Ø5.0 mm	5.0 mm	0.5 mm
45440K	Ø5.0 mm	5.0 mm	1.5 mm
45441K	Ø5.0 mm	5.0 mm	2.5 mm
45442K	Ø5.0 mm	5.0 mm	4.0 mm
45428K	Ø5.0 mm	6.0 mm	0.5 mm
45429K	Ø5.0 mm	6.0 mm	1.5 mm
45430K	Ø5.0 mm	6.0 mm	2.5 mm
45431K	Ø5.0 mm	6.0 mm	4.0 mm

**Quick-Abutment
Quick-Cap®**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
45444K	Ø5.0 mm, 1-Pack	5.0 mm
45434K	Ø5.0 mm, 1-Pack	6.0 mm

**Quick-Abutment
Temporary/Healing Cap**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
45237K	Ø5.0 mm	5.0 mm
45211K	Ø5.0 mm	6.0 mm

**Quick-Abutment
Analog**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
45449K	Ø5.0 mm	5.0 mm
45447K	Ø5.0 mm	6.0 mm



**Quick-Abutment
Plastic Sleeve**

Catalog #

Description

Flare

45241K
45243K
45217K
45220K

Ø5.0 mm, Locking
Ø5.0 mm, Non-Locking
Ø5.0 mm, Locking
Ø5.0 mm, Non-Locking

5.0 mm
5.0 mm
6.0 mm
6.0 mm



**Quick-Abutment
Ceramic Coping**

Catalog #

Description

Flare

45245K
45303K

Ø5.0 mm
Ø5.0 mm

5.0 mm
6.0 mm



**Reamer Tool for Quick-
Abutment Plastic Sleeve**

Catalog #

Description

45311K

Reamer Tool and Centering Pins




Prepable Abutment

Catalog #

Description

45044K

Ø5.0 mm




Ti-Base Abutment (Short)

Catalog #

Description

45808K

Ø5.0 mm



**Pink Ti-Base
Abutment (Short)**

Catalog #

Description

45009K

Ø5.0/5.5/6.5 mm



Ti-Base Abutment (Long)


Catalog #

Description

45903K

Ø5.0 mm

Multi Unit Abutment



Catalog

Description

Cuff Height

45319K	Ø5.0/5.5/6.5 mm, Straight	1.0 mm
45320K	Ø5.0/5.5/6.5 mm, Straight	2.0 mm
45939K	Ø5.0/5.5/6.5 mm, Straight	3.0 mm
45940K	Ø5.0/5.5/6.5 mm, Straight	4.0 mm
45941K*	Ø5.0/5.5/6.5 mm, Angled, 17°	3.0 mm
45942K*	Ø5.0/5.5/6.5 mm, Angled, 17°	4.0 mm

Multi Unit Abutment Healing Cap



Catalog

Description

45330K	Ø5.0/5.5/6.5 mm
---------------	-----------------

Multi Unit Abutment Impression Post



Catalog

Description

45327K	Open Tray, Ø5.0/5.5/6.5 mm with screw
45328K	Closed Tray, Ø5.0/5.5/6.5 mm

Multi Unit Abutment Analog



Catalog

Description

45329K	Ø5.0/5.5/6.5 mm
---------------	-----------------

Multi-Unit Abutment Ti Temporary Sleeve



Catalog

Description

45331K	Ø5.0/5.5/6.5 mm
---------------	-----------------

Multi Unit Abutment Plastic Sleeve



Catalog

Description

45332K	Ø5.0/5.5/6.5 mm
---------------	-----------------

Multi Unit Abutment Prosthetic Screw



Catalog

Description

45333K	1-Pack
---------------	--------

Replacement Abutment Screw for Multi Unit Angled Abutment




Catalog

Description

45334K	1-Pack
---------------	--------

*Available mid-2018



**LOCATOR R-Tx™
Abutment**



Catalog #

Description

Cuff Height

LOCR-TX 31200-01	Ø3.5/3.8 mm	1.0 mm
LOCR-TX 31200-02	Ø3.5/3.8 mm	2.0 mm
LOCR-TX 31200-03	Ø3.5/3.8 mm	3.0 mm
LOCR-TX 31200-04	Ø3.5/3.8 mm	4.0 mm
LOCR-TX 31200-05	Ø3.5/3.8 mm	5.0 mm
LOCR-TX 31200-06	Ø3.5/3.8 mm	6.0 mm
LOCR-TX 31201-01	Ø4.1/4.5 mm	1.0 mm
LOCR-TX 31201-02	Ø4.1/4.5 mm	2.0 mm
LOCR-TX 31201-03	Ø4.1/4.5 mm	3.0 mm
LOCR-TX 31201-04	Ø4.1/4.5 mm	4.0 mm
LOCR-TX 31201-05	Ø4.1/4.5 mm	5.0 mm
LOCR-TX 31201-06	Ø4.1/4.5 mm	6.0 mm
LOCR-TX 31202-01	Ø5.0/5.5/6.5 mm	1.0 mm
LOCR-TX 31202-02	Ø5.0/5.5/6.5 mm	2.0 mm
LOCR-TX 31202-03	Ø5.0/5.5/6.5 mm	3.0 mm
LOCR-TX 31202-04	Ø5.0/5.5/6.5 mm	4.0 mm
LOCR-TX 31202-05	Ø5.0/5.5/6.5 mm	5.0 mm
LOCR-TX 31202-06	Ø5.0/5.5/6.5 mm	6.0 mm

**LOCATOR R-Tx™
Abutment Accessories**



Catalog #

Description

LOCR-TX 30013-01	Denture Attachment Processing Assembly, 4-Pack
LOCR-TX 30001-01	Zero Retention Insert, 4-Pack
LOCR-TX 30002-01	Low Retention Insert, 4-Pack
LOCR-TX 30003-01	Medium Retention Insert, 4-Pack
LOCR-TX 30004-01	High Retention Insert, 4-Pack
LOCR-TX 30005-01	Extra-High Retention Insert, 4-Pack
LOCR-TX 30014-01	3.35 mm Abutment Analog, 4-Pack
LOCR-TX 30015-01	4.0 mm Abutment Analog, 4-Pack
LOCR-TX 30016-01	5.0 mm Abutment Analog, 4-Pack
LOCR-TX 30017-01	Impression Coping, 4-Pack
LOCR-TX 30021-01	Retention Insert Tool
LOCR-TX 04936	.050"/1.25 mm Hex, Torque Wrench
LOCR-TX 08007	.050"/1.25 mm Hex Driver with Thumb Knob
LOCR-TX 08008	.050"/1.25 mm Hex, Latch Lock Driver
LOCR-TX 30012-01	Processing Insert, 4-Pack
LOCR-TX 30018-01	Processing Spacer, 4-Pack




LOCATOR® Abutment Assembly

<i>Catalog #</i>	<i>Description</i>	<i>Cuff Height</i>
45062K	Ø3.5 mm	1.0 mm
45258K	Ø3.5 mm	2.0 mm
45259K	Ø3.5 mm	3.0 mm
45260K	Ø3.5 mm	4.0 mm
45063K	Ø3.5 mm	5.0 mm
45064K	Ø3.5 mm	6.0 mm
45065K	Ø4.1 mm	1.0 mm
45261K	Ø4.1 mm	2.0 mm
45262K	Ø4.1 mm	3.0 mm
45263K	Ø4.1 mm	4.0 mm
45066K	Ø4.1 mm	5.0 mm
45067K	Ø4.1 mm	6.0 mm
45068K	Ø5.0 mm	1.0 mm
45264K	Ø5.0 mm	2.0 mm
45265K	Ø5.0 mm	3.0 mm
45069K	Ø5.0 mm	4.0 mm
45070K	Ø5.0 mm	5.0 mm
45071K	Ø5.0 mm	6.0 mm

LOCATOR® Abutment Accessories



<i>Catalog #</i>	<i>Description</i>
L85192	Denture Cap with Black Processing Male, 2-Pack
45266K	LOCATOR® Replacement Male, White, Standard Retention, 4-Pack
45267K	LOCATOR® Replacement Male, Pink, Light Retention, 4-Pack
45268K	LOCATOR® Replacement Male, Blue, (ELR) Extra Light Retention, 4-Pack
45269K	LOCATOR® Replacement Male, Green, (ER) Extended Range, 4-Pack
45270K	LOCATOR® Female Analog, 4-Pack
45271K	LOCATOR® Impression Coping, 4-Pack
45273K	LOCATOR® Core Tool
45274K	LOCATOR® Torque Wrench Driver



Catalog # Surgical Options

15278K	Surgical Cassette (empty)
15241K	A. PrimaConnex® Surgical Kit
15244K	B. PrimaConnex® Tapered Drill and Tap Set
15249K	C. PrimaConnex® Straight Drill and Tap Set

The PrimaConnex® Surgical Kit contains all instrumentation except the final tapered and/or straight drills and taps. The Tapered and/or Straight Drill and Tap Set(s) must be ordered separately, depending on the implant design to be placed.

Catalog # Description

15241K	A. PrimaConnex® Surgical Kit <i>Includes:</i>
15278K	Surgical Cassette
15200K	1.5 x 15 mm Initial Drill
15201K	1.8 x 11.5 mm Initial Drill
15202K	1.8 x 15 mm Initial Drill
15225K	Depth Probe
L1404-01K	Drill Extender
15286K	3.5/3.8 Implant Driver, TiLobe Connection, Short
15287K	3.5/3.8 Implant Driver, TiLobe Connection, Long
15288K	4.1/4.5 Implant Driver, TiLobe Connection, Short
15289K	4.1/4.5 Implant Driver, TiLobe Connection, Long
15290K	5.0/5.5/6.5/7.0/8.0/9.0 Implant Driver, TiLobe Connection, Short
15282K	3.5 Implant Driver, PrimaConnex®, Ratchet
15283K	4.1 Implant Driver, PrimaConnex®, Ratchet
15284K	5.0 Implant Driver, PrimaConnex®, Ratchet

Catalog # Description

L1800-22-02K	Irrigation Port Cleaner
15226K	Parallel Pins, Straight, 1.5 mm/1.8 mm, 2-Pack
45045K	Parallel Pins, 15° Angled, 1.5 mm/1.8 mm, 2-Pack
15274K	Quad Driver, Latch, Short
15275K	Quad Driver, Latch, Long
15272K	Quad Driver, Swivel, Short
15273K	Quad Driver, Swivel, Long
15224K	Ratchet Adapter
L1702-01K	Surgical Ratchet
L1407K	Conversion Handle, Latch
4409-0255	PrimaConnex® Surgical Manual
15240K	Implant Sizing Overlays
4409-0257	Wall Chart
4409-0258	Drill Usage Chart



Catalog #	Description
15244K	B. PrimaConnex® Tapered Drill and Tap Set <i>Includes:</i>
15205K	3.5 x 10 mm Final Tapered Drill
15206K	3.5 x 11.5 mm Final Tapered Drill
15207K	3.5 x 13 mm Final Tapered Drill
15208K	3.5 x 15 mm Final Tapered Drill
15209K	4.1 x 10 mm Final Tapered Drill
15210K	4.1 x 11.5 mm Final Tapered Drill
15211K	4.1 x 13 mm Final Tapered Drill
15212K	4.1 x 15 mm Final Tapered Drill
15213K	5.0 x 10 mm Final Tapered Drill
15214K	5.0 x 11.5 mm Final Tapered Drill
15215K	5.0 x 13 mm Final Tapered Drill
15075K	5.0 x 15 mm Final Tapered Drill
15217K	3.5 x 10/11.5 mm Tapered Surgical Tap
15247K	3.5 x 13/15 mm Tapered Surgical Tap
15218K	4.1 x 10/11.5 mm Tapered Surgical Tap
15248K	4.1 x 13/15 mm Tapered Surgical Tap
15219K	5.0 x 10/11.5 mm Tapered Surgical Tap
15076K	5.0 x 13/15 mm Tapered Surgical Tap

Catalog #	Description
15249K	C. PrimaConnex® Straight Drill and Tap Set <i>Includes:</i>
15250K	2.8 mm x 11.5 mm Straight Drill
15251K	2.8 mm x 15 mm Straight Drill
15252K	3.3 mm x 11.5 mm Straight Drill
15253K	3.3 mm x 15 mm Straight Drill
15254K	4.2 mm x 11.5 mm Straight Drill
15255K	4.2 mm x 13 mm Straight Drill
15101K	1.8 mm/2.8 mm Step Drill
15256K	2.8 mm/3.3 mm Step Drill
15257K	3.3 mm/4.2 mm Step Drill
15258K	3.3 mm Countersink Drill
15259K	4.0 mm Countersink Drill
15260K	5.0 mm Countersink Drill
15262K	3.3 x 15 mm Straight Surgical Tap
15264K	4.0 x 15 mm Straight Surgical Tap
15265K	5.0 x 13 mm Straight Surgical Tap



Tissue Punch

Catalog

Description

15237K	Ø3.5 mm, with Guide Pin
15238K	Ø4.1 mm, with Guide Pin
15239K	Ø5.0 mm, with Guide Pin




Round Marking Bur

Catalog

Description

L1405-23-00K	Ø2.3 mm
---------------------	---------




Initial Drills

Catalog

Description

15200K	Ø1.5 x 15 mm, Non-Irrigated
15201K	Ø1.8 x 11.5 mm, Non-Irrigated
15202K	Ø1.8 x 15 mm, Non-Irrigated



Parallel Pins

Catalog

Description

15226K
45045K

Ø1.5/1.8 mm, Straight, 2-Pack
Ø1.5/1.8 mm, 15°, 2-Pack



Final Tapered Drills

Catalog

Description

15203K
15204K

Ø3.0 x 13 mm
Ø3.0 x 15 mm

15205K
15206K
15207K
15208K

Ø3.5 x 10 mm
Ø3.5 x 11.5 mm
Ø3.5 x 13 mm
Ø3.5 x 15 mm

15209K
15210K
15211K
15212K

Ø4.1 x 10 mm
Ø4.1 x 11.5 mm
Ø4.1 x 13 mm
Ø4.1 x 15 mm

15213K
15214K
15215K
15075K

Ø5.0 x 10 mm
Ø5.0 x 11.5 mm
Ø5.0 x 13 mm
Ø5.0 x 15 mm



Removable Drill Stop for Tapered Drills

Catalog

Description

15220K
15221K
15222K
15223K

3.0 mm
3.5 mm
4.1 mm
5.0 mm



Tapered Surgical Taps

Catalog

Description

15217K
15247K

Ø3.5 x 10/11.5 mm
Ø3.5 x 13/15 mm

15218K
15248K

Ø4.1 x 10/11.5 mm
Ø4.1 x 13/15 mm

15219K
15076K

Ø5.0 x 10/11.5 mm
Ø5.0 x 13/15 mm

Only PrimaConnex® Initial Drills are non-irrigated. All other Prima Tapered Drills, Tapered Taps, Straight Drills, Step Drills, Countersink Drills and Straight Taps are irrigated.



Straight Drills



Catalog

Description

15250K	Ø2.8 x 11.5 mm
15251K	Ø2.8 x 15 mm
15252K	Ø3.3 x 11.5 mm
15253K	Ø3.3 x 15 mm
15254K	Ø4.2 x 11.5 mm
15255K	Ø4.2 x 13 mm

Step Drills for Straight Implants



Catalog

Description

15101K	Ø1.8/Ø2.8 mm
15256K	Ø2.8/Ø3.3 mm
15257K	Ø3.3/Ø4.2 mm

Countersink Drills for Straight Implants



Catalog

Description

15258K	For Ø3.3 mm Straight Implants
15259K	For Ø4.0 mm Straight Implants
15260K	For Ø5.0 mm Straight Implants

Straight Surgical Taps



Catalog

Description

15261K	Ø3.3 x 11.5 mm
15262K	Ø3.3 x 15 mm
15263K	Ø4.0 x 11.5 mm
15264K	Ø4.0 x 15 mm
15265K	Ø5.0 x 13 mm

Drill Extender



Catalog

L1404-01K

Tissue Punch



Catalog

Description

15236K	Ø3.0 mm Tissue Punch with Guide Pin
15237K	Ø3.5 mm Tissue Punch with Guide Pin
15238K	Ø4.1 mm Tissue Punch with Guide Pin
15239K	Ø5.0 mm Tissue Punch with Guide Pin

Bone Profiler



Catalog


Description

15279K	Ø3.5 mm, with Guide Pin
15280K	Ø4.1 mm, with Guide Pin
15281K	Ø5.0 mm, with Guide Pin

Only PrimaConnex® Initial Drills are non-irrigated. All other Prima Tapered Drills, Tapered Taps, Straight Drills, Step Drills, Countersink Drills and Straight Taps are irrigated.



**PrimaConnex®
Implant Driver, Ratchet**



Catalog #

Description

15282K
15283K
15284K

Ø3.5 mm
Ø4.1 mm
Ø5.0 mm

**PrimaConnex®
Implant Driver, Latch**



Catalog #

Description

Length

15286K
15287K
15288K
15289K
15290K

Ø3.5/3.8 mm Implant Driver, Short
Ø3.5/3.8 mm Implant Driver, Long
Ø4.1/4.5 mm Implant Driver, Short
Ø4.1/4.5 mm Implant Driver, Long
Ø5.0/5.5/6.5/7.0/8.0/9.0 mm Implant Driver, Short

16 mm
21 mm
16 mm
21 mm
16 mm

Surgical Ratchet



Catalog #

Description

L1702-01K

Surgical Ratchet

**Surgical Ratchet/
Torque Wrench**



Catalog #

Description

K15695

Surgical Ratchet/Torque Wrench

Quad Drivers



Catalog #

Description

15272K
15273K

Swivel Head, Short
Swivel Head, Long

15274K
15275K

Latch Type, Short
Latch Type, Long

15276K
15277K

Torque Wrench, Short
Torque Wrench, Long



Ratchet Adapter

Catalog #

15224K



Directional Indicators

Catalog #

45041K

Description

Try-In Set




Conversion Handles

Catalog #

L1407K

Description

Conversion Handle, Latch Type



Surgical Hand Driver

Catalog #

15271K



Depth Probe

Catalog #

15225K




Irrigation Port Cleaner

Catalog #

L1800-22-02K



Surgical Kit Overview



Catalog #	Description
15651K	Tapered PrimaConnex® Guided Surgery Kit (with instruments) Includes the instrumentation required to place PrimaConnex® Tapered Implants
15650K	Guided Surgery Cassette (without instruments)

Anchor Pins & Anchor Pin Drill



Catalog

E40100
E20101

Description

Anchor Pins (package of 3)
Anchor Pin Drill Ø1.6 mm

Tissue Punch



Catalog

15671K
15672K
15673K

Description

PrimaConnex® Guided Tissue Punch, Ø3.5
PrimaConnex® Guided Tissue Punch, Ø4.1
PrimaConnex® Guided Tissue Punch, Ø5.0

Guided Surgery Drills



Catalog

15674K
15675K

Description

PrimaConnex® Guided Initial Drill, Short
PrimaConnex® Guided Initial Drill, Long

15652K
15653K
15654K
15655K

PrimaConnex® Guided Surgery, Tapered Drill, Ø3.5 x 10
PrimaConnex® Guided Surgery, Tapered Drill, Ø3.5 x 11.5
PrimaConnex® Guided Surgery, Tapered Drill, Ø3.5 x 13
PrimaConnex® Guided Surgery, Tapered Drill, Ø3.5 x 15

15656K
15657K
15658K
15659K


PrimaConnex® Guided Surgery, Tapered Drill, Ø4.1 x 10
PrimaConnex® Guided Surgery, Tapered Drill, Ø4.1 x 11.5
PrimaConnex® Guided Surgery, Tapered Drill, Ø4.1 x 13
PrimaConnex® Guided Surgery, Tapered Drill, Ø4.1 x 15

15660K
15661K
15662K
15663K

PrimaConnex® Guided Surgery, Tapered Drill, Ø5.0 x 10
PrimaConnex® Guided Surgery, Tapered Drill, Ø5.0 x 11.5
PrimaConnex® Guided Surgery, Tapered Drill, Ø5.0 x 13
PrimaConnex® Guided Surgery, Tapered Drill, Ø5.0 x 15



Drill Guides



Catalog

15676K
15678K
15679K


Description

PrimaConnex® Guided Drill Guide, Initial Drill, Ø3.5
PrimaConnex® Guided Drill Guide, Initial Drill, Ø4.1
PrimaConnex® Guided Drill Guide, Initial Drill, Ø5.0

15680K
15681K
15682K

PrimaConnex® Guided Drill Guide, Ø3.5/Ø4.1
PrimaConnex® Guided Drill Guide, Ø3.5/Ø5.0
PrimaConnex® Guided Drill Guide, Ø4.1/Ø5.0

Guided Surgery Bone Taps




Catalog

15664K
15665K
15667K
15668K
15669K
15670K

Description

PrimaConnex® Guided Bone Tap, Ø3.5, Short
PrimaConnex® Guided Bone Tap, Ø3.5, Long
PrimaConnex® Guided Bone Tap, Ø4.1, Short
PrimaConnex® Guided Bone Tap, Ø4.1, Long
PrimaConnex® Guided Bone Tap, Ø5.0, Short
PrimaConnex® Guided Bone Tap, Ø5.0, Long

Guided Implant Drivers



Catalog

15683K
15684K
15685K

Description

PrimaConnex® Guided Implant Driver, Ø3.5
PrimaConnex® Guided Implant Driver, Ø4.1
PrimaConnex® Guided Implant Driver, Ø5.0

Surgical Ratchet & Ratchet Adapter




Catalog

K15695
15686K

Description

Surgical Ratchet/Torque Wrench
Ratchet Adapter

Quad Drivers



Catalog

15272K
15273K
15276K
15277K

Description

Quad Swivel, Short
Quad Swivel, Long
Quad Driver, Torque Wrench, Short
Quad Driver, Torque Wrench, Long

Terms and Conditions

TERMS AND CONDITIONS OF SALE. The following are the terms and conditions by which Keystone Dental, Inc. sells its products. As used herein, the term “Keystone” means Keystone Dental, Inc. and its applicable affiliates; the term “Products” means the Keystone products and services described in this invoice, and “Purchaser” means the person or company purchasing the Products.

GOVERNING TERMS. Any purchase and shipment of Products shall be deemed to be on the terms and conditions stated herein. Any and all terms and conditions submitted by Purchaser are hereby rejected.

TAXES. The prices set forth herein do not include any sales, use, excise, ad valorem, property or other taxes applicable to the sale, use or delivery of the Products, all of which shall be paid by Purchaser separately or added to the contract price and paid by Purchaser to Keystone.

PRICES. The prices set forth for the Products are in United States Dollars. All prices shall be as of the invoice date, regardless of the date the Product is originally ordered from Keystone.

DELIVERY. Products sold hereunder shall be delivered F.O.B. shipping point. Purchaser shall assume responsibility for all subsequent delivery or shipping charges.

ACCEPTANCE. Purchaser must inspect all Products for defects and/or non-conformity within seven (7) calendar days of delivery. After this period, Purchaser will have accepted such Products.

PAYMENT. Payment terms shall be as set forth on the front of the invoice. If payment is late, the maximum statutory interest rate will apply on such late amount. Keystone may, in its sole discretion, suspend shipments of Products until full payment by Purchaser of all past due amounts to Keystone is made.

RETURN; EXCHANGE. Purchaser may return any Product within thirty (30) calendar days of the invoice date so long as there is an offsetting (exchange) order of equal or greater value. Purchaser may return any Product after thirty (30) calendar days from the invoice date, but no later than ninety (90) calendar days after the invoice date, so long as there is an offsetting (exchange) order of equal or greater value, provided, however, such returns shall be subject to a 15% re-stocking fee. Keystone will not accept any returns or exchanges after ninety (90) calendar days following the invoice date.

No biomaterial products (bone grafting and membrane) will be accepted for return or exchange. No returns or exchanges will be accepted for Products that (i) were originally sold as part of a kit, pack or set, (ii) were included on the original invoice at no charge, (iii) have less than six (6) months remaining before expiration or (iv) were custom ordered. The Purchaser shall pay all shipping costs associated with the return or exchange of a Product. No Products will be accepted by Keystone for return or exchange without an RMA number given to Purchaser by Keystone Customer Service. Returned Products must be unexpired, unopened and undamaged, returned freight prepaid.

WARRANTY. EXCEPT AS EXPRESSLY SET FORTH IN KEYSTONE'S LIFETIME IMPLANT WARRANTY PROGRAM, KEYSTONE DOES NOT MAKE ANY WARRANTY, EXPRESS OR IMPLIED, EXCEPT THAT ITS PRODUCTS SHALL BE FREE FROM DEFECTS IN MATERIAL AND/OR WORKMANSHIP. THIS WARRANTY APPLIES ONLY TO THE ORIGINAL PURCHASER. PURCHASER ASSUMES ALL RISKS AND LIABILITY RESULTING FROM THE USE OF THESE PRODUCTS, WHETHER USED SEPARATELY OR IN COMBINATION WITH OTHER PRODUCTS. THE WARRANTY DOES NOT APPLY IF THE FAILURE OF THE PRODUCT IS CAUSED BY TRAUMA OR ACCIDENT OR BY THE PATIENT IN WHOM THE IMPLANT WAS INSTALLED OR IF THE PURCHASER IS NOT CURRENT ON ALL OUTSTANDING PAYMENTS TO KEYSTONE AT THE TIME OF SUCH CLAIM. KEYSTONE MAY MODIFY OR TERMINATE ITS WARRANTY PROGRAM AT ANY TIME, HOWEVER ANY CHANGES TO THE WARRANTY PROGRAM WILL NOT AFFECT PRODUCTS ALREADY PURCHASED AT THE TIME OF SUCH CHANGE. PURCHASER'S SOLE REMEDY WITH RESPECT TO CLAIMS RELATED TO KEYSTONE PRODUCTS IS RESTRICTED TO THE WARRANTY POLICY.

LIMITATION OF LIABILITY. UNDER NO CIRCUMSTANCES SHALL KEYSTONE OR ANY OF ITS AFFILIATES BE LIABLE FOR BUSINESS INTERRUPTION, LOSS OF PROFITS, SPECIAL, INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE AND FROM ANY CAUSE WHETHER BASED IN CONTRACT, WARRANTY OR TORT (INCLUDING NEGLIGENCE) OR OTHER LEGAL THEORY, EVEN IF KEYSTONE HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

INTELLECTUAL PROPERTY. All intellectual property rights in Products sold to Purchaser shall remain the exclusive property of Keystone and its licensors.

PERSONAL DATA. It is Purchaser's responsibility to obtain patient's consent to transfer any personal data to Keystone.

TERMINATION. Keystone may terminate this or any agreement with Purchaser to purchase its products upon immediate written notice if Purchaser fails to pay any of its debt obligations to Keystone on time.

FORCE MAJEURE. Keystone shall not be liable for any loss or damage due to failure or delay arising out of any cause beyond the reasonable control of Keystone. In the event of any failure or delay resulting from such causes, an equitable adjustment of delivery and any other appropriate terms and conditions shall be made. No such failure or delay shall be the basis for an increase in Keystone's obligations, nor any termination by Purchaser.

VALIDITY. If any provision of these Terms and Conditions is found to be illegal or unenforceable in any respect, such illegality or unenforceability shall not affect any other provision of these Terms and Conditions, all of which shall remain enforceable in accordance with their terms.

GENERAL PROVISIONS. The purchase of the Products and these Terms and Conditions are governed by the laws of the Commonwealth of Massachusetts, without regard to conflicts of law principles. Any dispute, claim or controversy arising out of or in connection with the purchase of the Products shall be resolved by arbitration in Boston, Massachusetts. This document incorporates all oral and written representations between the parties and constitutes the entire agreement and understanding of the parties with respect to the subject matter hereof and supersedes any and all other agreements either oral or written between the parties with respect to such subject matter. No amendment or modification may be made to this document unless in writing and duly executed by an authorized representative of Keystone. Keystone may subcontract, assign or transfer its obligations or rights to a competent third party. Purchaser may not assign or transfer any of its obligations without the consent of Keystone.

Keystone Dental, Inc.
Global Headquarters

154 Middlesex Turnpike
Burlington, MA 01803 USA
General Inquiries: 781-328-3600
Customer Relations: 866-902-9272 or 781-328-3490
Fax: 866-903-9272 or 781-328-3400
Email: info@keystonedental.com



Caution, consult accompanying documents

Trademark Acknowledgements

TiLobe, PrimaPlus, Quick-Cap and "Smarter Thinking. Simpler Design." are trademarks and registered trademarks of Keystone Dental, Inc.

Vertrieb / Distribution

Deutschland, Belgien,
Niederlande, Luxemburg,
Frankreich:

RUNDAS GmbH

Amalienstr. 62
46537 Dinslaken

Tel.: 02064 625 95 50
Fax: 02064 625 95 80

E-Mail: info@rundas.de
Internet: www.rundas.de